Nea Medical Device

04/2020

47206 / 2

Designed in the Netherlands

www.push.eu

The Netherlands

6199 AB Maastricht-Airport

Nea International bv


4.30.2

REF

# PUSH SPORTS

## PATELLA BRACE

### Useful information for you!

## ENGLISH

### PUSH SPORTS PATELLA BRACE

The viscoelastic pad in the Push Sports Patella Brace applies pressure to the patellar tendon immediately below the kneecap and reduces pain. Pressure on this tendon provides relief for the cartilage on the underside of the kneecap, and the tendon insertion point at the front. The brace is easy to apply, and to adjust to the tension that you require. This enables you to set the degree of compression, while the brace stays securely in place, partly because of the anti-slip properties of the pad.

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- INDICATIONS**
- Tendinitis of the patellar ligament
  - Jumper's knee
  - Patellofemoral pain syndrome
  - Osgood-Schlatter disease

- CONTRAINDICATIONS**
- There are no known disorders where the Push Sports Patella Brace should not be worn.

- CAUTIONS**
- Ensure that the brace is not fitted too tightly. Use the brace only for the appropriate indications. Proper use of your Push Sports Brace may require diagnosis by a doctor or therapist. Pain is always a warning sign. If pain or complaints persist without a clear cause, it is recommended that you visit your GP, physiotherapist or sports physician.

- Push Sports braces were designed to be worn on an undamaged skin. Damaged skin must first be covered with suitable material.
- Optimal functioning can only be achieved by selecting the right size. Therefore, preferably fit the brace before purchasing, to check its size. In the case of an online purchase, always ensure that there is a good return service. Your reseller's professional staff will explain the best way to fit your Push Sports brace. We recommend that you keep this product information in a safe place, to be able to consult it for fitting instructions whenever needed.
- Check your Push Sports brace each time before use for signs of wear or ageing of components or seams. Optimal performance of your Push Sports brace can only be guaranteed if the product is fully intact.

- Within the framework of EU regulation 2017/745 on medical devices, the user must report any serious incident involving the product to the manufacturer and to the competent authority of the Member State in which the user is established.

### SPECIAL FEATURES

The patella brace pushes against the patella tendon immediately below the patella, causing a slight tilting and relieving the cartilage at the back of the patella. During each step, the brace provides supportive pressure on the patella tendon, absorbing shock and reducing the impact on the patella.

### TAKING MEASUREMENTS

The patella brace is available in one size, universally applicable for the left and right sides.

### FITTING

**Before use:** remove the protective foil from the pad.

- 1 Bend the leg. Place the patella brace on the knee tendon immediately below the kneecap, with the kneecap fitting in the recess.
- 2 Hold the pad with one hand and use the other hand to wrap the brace around the leg.
- 3 Close the hook and loop material as tightly as possible, without pinching. If the strap fails to close on the hook and loop material, use the strap's other fastener.

[pushsports.eu/patella-brace/applying](https://pushsports.eu/patella-brace/applying)

For details of the distributor in your country, we refer to [www.push.eu/distributors](http://www.push.eu/distributors)

### WASHING INSTRUCTIONS/MAINTENANCE

To promote hygiene and to maximise the lifespan of the brace, it is important to wash the brace regularly, observing the washing instructions. The brace can be washed at temperatures of up to thirty degrees Celsius, on a delicate fabric cycle or by hand. Before washing, fasten the hook and loop material to protect the brace and other laundry. Do not use any bleaching agents. Let the brace dry in the open air (not near a source of heat or in a tumble drier). Never store the brace when it is wet. Always let it dry first.

- IMPORTANT**
- The Push Sports Patella Brace is not a replacement for knee protectors used in indoor sports.

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