Material composition / Materialzusammensetzung	
Polyamide / Polyamid	44%
Polyester / Polyester	31%
Polyurethane / Polyurethan	16%
Elasthane / Elastan	6%
Others / Sonstige	3%

SUPORTE PARA O PULSO PUSH SPORTS

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PUSH SPORTS WRIST SUPPORT

The Push Sports Wrist Support provides compressive support around the wrist. This reduces pain, and provides effective support and a sense of security to control minor strains during sports and exercise. The Push Sports Wrist Support is manufactured from high-quality materials that ensure comfort and adjustability to each user's specific needs.

WHEN TO USE THE PUSH SPORTS WRIST SUPPORT?

The Push Sports Wrist Support offers relief in the event of wrist irritations involving minor pain or a sense of mild instability. Use of the wrist support ensures ligaments and/or tendons in the wrist joint are protected against strain during exercise and sports, such as gymnastics, yoga or tennis.

INDICATIONS

- To prevent recurring strain symptoms in the wrist joint
- Mild sprains in the wrist joint
- Laxity in the wrist joint

CONTRAINDICATIONS

- There are no known disorders where the Push Sports Wrist Support should not be worn.

CAUTIONS

- Ensure that the brace is not fitted too tightly. Use the brace only for the appropriate indications. Proper use of your Push Sports brace may require diagnosis by a doctor or therapist. Pain is always a warning sign. If pain or complaints persist without a clear cause, it is recommended that you visit your GP, physiotherapist or sports physician.
- Push Sports braces were designed to be worn on an undamaged skin. Damaged skin must first be covered with suitable material.
- Optimal functioning can only be achieved by selecting the right size. Therefore, preferably fit the brace before purchasing, to check its size. In the case of an online purchase, always ensure that there is a good return service. Your reseller's professional staff will explain the best way to fit your Push Sports brace. We recommend that you keep this product information in a safe place, to be able to consult it for fitting instructions whenever needed.

- Check your Push Sports brace each time before use for signs of wear or ageing of components or seams. Optimal performance of your Push Sports brace can only be guaranteed if the product is fully intact.
- Within the framework of EU regulation 2017/745 on medical devices, the user must report any serious incident involving the product to the manufacturer and to the competent authority of the Member State in which the user is established.

SPECIAL FEATURES

The Push Sports Wrist Support provides ulnar compression that relieves pain in the wrist. The brace supports the wrist with very slight restriction of movement.

TAKING MEASUREMENTS

The wrist support is available in one size, with left and right versions.

FITTING

- Insert the thumb in the thumb loop, with the Push Sports logo visible.
- Pull the support across the palm to the outside of the hand, across the back of the hand, along the front again, and fasten the hook and loop material.
- Remove the loop from the thumb and fasten the strap on the hook and loop material at the back of the hand (Fig. ④).

pushsports.eu/wrist-support/applying

For details of the distributor in your country, we refer to www.push.eu/distributors

WASHING INSTRUCTIONS/MAINTENANCE

To promote hygiene and to maximise the lifespan of the brace, it is important to wash the brace regularly, observing the washing instructions. The brace can be washed at temperatures of up to thirty degrees Celsius, on a delicate fabric cycle or by hand. Before washing, fasten the hook and loop material to protect the brace and other laundry. Do not use any bleaching agents. Let the brace dry in the open air (not near a source of heat or in a tumble drier). Never store the brace when it is wet. Always let it dry first.

WRIST SUPPORT

Useful information for you!

Nea International bv
Europalaan 31
6199 AB Maastricht-Alroort
The Netherlands
www.push.eu
Made in the EU
47201 / 2
04/2020

Medical Device
MD
REF

Push Sports
Wrist Support

Push Sports is a registered trademark of Nea International bv